

# SEQUENCE LISTING

- <110> Hyeyoung Lee, Hye Eun Bang, Sang-Nae Cho, Gill-Han BAI,  
Sang-Jae Kim
- <120> A method for identifying Micobacteria tuberculosis and  
non-tuberculosis Micobacteria, together with detecting resistance  
to an antituberculosis drug of Micobacteria obtained by mutation  
of rpoB gene
- <130> 0217-0008
- <160> 30
- <170> KopatentIn 1.71
- <210> 1
- <211> 25
- <212> DNA
- <213> Artificial Sequence
- <220>
- <223> MOTT-rpo-long-B-5' primer for PCR amplication of rpoB gene
- <400> 1  
tcaaggagaa gcgctacgac ctggc 25
- <210> 2
- <211> 24
- <212> DNA
- <213> Artificial Sequence
- <220>
- <223> TR8-long-NB-3' primer for PCR amplication of rpoB gene
- <400> 2  
acgggtgcac gtcgcggacc tcca 24
- <210> 3
- <211> 20
- <212> DNA
- <213> Artificial Sequence
- <220>
- <223> Oligomer probe for all types of Mycobacteria
- <400> 3  
gacgtcgtcg ccaccatcga 20
- <210> 4

<211> 15  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Oligomer probe for M. tuberculosis complex

<400> 4  
 catgtcggcg agccc 15

<210> 5  
 <211> 19  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Oligomer probe for M. avium

<400> 5  
 aaacggtgag ccgatcacc 19

<210> 6  
 <211> 18  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Oligomer probe for M. intracellulae

<400> 6  
 aaacctgcac gcgggcga 18

<210> 7  
 <211> 21  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Oligomer probe for M. scrofulaceum

<400> 7  
 aaaaacgtac ggatggccag c 21

<210> 8  
 <211> 19  
 <212> DNA  
 <213> Artificial Sequence

<220>



<210> 13  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Oligomer probe for *M. chelonae*

<400> 13  
 aaatggtgac tgccaccacg 20

<210> 14  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Oligomer probe for *M. abscessus*

<400> 14  
 aaaaggtgac caccaccacc 20

<210> 15  
 <211> 15  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Oligomer probe for *M. ulcerans*

<400> 15  
 ggccagccca tcacc 15

<210> 16  
 <211> 16  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Oligomer probe for *M. genavense*/*M. simiae*

<400> 16  
 ccagccgacg atgacg 16

<210> 17  
 <211> 19  
 <212> DNA

<213> Artificial Sequence

<220>

<223> Oligomer probe for *M. gordonae* type I, III, IV

<400> 17

aaagtcggcg atccgatca

19

<210> 18

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> Oligomer probe for *M. gordonae* type II

<400> 18

aaaaacgtcg gcaagccga

19

<210> 19

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> Oligomer probe for *M. szulgai*

<400> 19

aaatctgaac gtcggcgag

19

<210> 20

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> Oligomer probe for *M. terrae*

<400> 20

aaagctcagg acggtcagt

19

<210> 21

<211> 18

<212> DNA

<213> Artificial Sequence

<220>

<223> Oligomer probe for Wild Type 509-514

<400> 21  
aaccagctga gccaatc 18

<210> 22  
<211> 18  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for M. Wild Type 515-520

<400> 22  
atggaccaga acaaccg 18

<210> 23  
<211> 18  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for Wild Type 521-525

<400> 23  
aaactgtcgg ggttgacc 18

<210> 24  
<211> 18  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for Wild Type 524-529

<400> 24  
ttgaccacaca agcgccga 18

<210> 25  
<211> 16  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for Wild Type 530-534

<400> 25  
ctgtcggcgc tggggc 16

<210> 26  
<211> 16  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for Mutant Type 531TTG

<400> 26  
ctgttggcgc tggggc 16

<210> 27  
<211> 18  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for Mutant Type 526 AAC

<400> 27  
aaaaccaaca agcgccga 18

<210> 28  
<211> 19  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for Mutant Type 516 GTC

<400> 28  
aatggtccag aacaaccg 19

<210> 29  
<211> 18  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for Mutant Type 513 CCA

<400> 29  
aaagctgacc ccattcat 18

<210> 30  
<211> 18  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> Oligomer probe for Mutant Type 511CCG

<400> 30  
aaagccgagc ccattcat

18

2005-10-23 10:50:10